
onsen

10
YEAR
WARRANTY

User Guide



Gebrauchs-anweisung | Instrucciones del usuario
Mode d'emploi | Istruzioni per l'uso
Gebruikershandleiding | Instrukcja obsługi

A NEW GENERATION
of open-cell foam
by **ONSEN**

ENSO

Poduszka ortopedyczna
Orthopedic Pillow
Orthopädisches Kissen

Born to sleep

MEMORY FOAM**SINKING OF THE HEAD****NO PROPER SUPPORT
FOR THE CERVICAL SPINE****CLOSED-CELL FOAM BLOCKS AIRFLOW,
CAUSING OVERHEATING****ONSEN FOAM****HEAD IS PROPERLY SUPPORTED****CORRECT POSITION
OF THE CERVICAL SPINE****OPEN-CELL FOAM
ALLOWS FREE AIRFLOW****Indications:**

- Arthritic pain of the cervical spine, discopathic changes of the spine
- Cervical hyperlordosis
- Non-specific nape and neck pain, certain types of headaches
- Degeneration of the intervertebral joints of the cervical spine

Advantages:

It causes no overheating and sweating of the head, which was the main disadvantage of visco-elastic foams used in pillows so far (the so-called memory foams/thermoelastic foams).

Pillow properties:

Owing to its specially designed shape, the pillow adapts to the curvature of the neck section of the spine, keeping it in a neutral, physiological position. This has a positive effect on the prevention of various diseases of the cervical spine. The core of the pillow is made of **ONSEN** high resilience foam with open cells. With these features it both enables excellent stabilisation and ensures free airflow, which guarantees a good night's sleep.

This is a **medical device** intended for use by non-professional users. You must use it according to this user manual. The product should be used by one patient. If you encounter any difficulties with fitting or using the product, please contact the manufacturer or an orthopedic doctor. The cushion has no contraindications for use. Its use does not require training. In case of worsening symptoms, discontinue use and immediately contact an orthopedic doctor. The product is registered in the EUDAMED database (PL-MF-000037523).

Pillow cleaning:**Foam:**

The pillow core (foam) should be hand washed at maximum 30 °C using a mild detergent. After washing, the accumulated water should be gently squeezed out (do not wring, fold or crease the pillow) and then the pillow should be left to dry at room temperature.

The foam should not be washed in a washing machine nor dried in a dryer or washer-dryer.

Pillowcase:

The pillowcase may be washed in a washing machine at maximum 40 °C using standard washing agents.

Storage:

The pillow must not be folded or creased and must be stored and transported unfolded. Store it in a dry place.

Do not expose the pillow to direct sunlight.

Do not expose the pillow to high or low temperatures (above 40 °C and below 8 °C).


30°C

40°C

UNIVERSAL SIZE

Pillow dimensions:

width 62 cm
 length 38 cm
 max. height 11 cm
 (tolerance +/- 2 cm)

Sleeping positions on a double-sided pillow ENSO


Back position


When sleeping on your back, your head should rest in the concave section in the centre of the pillow. In this position it is raised 3 to 4 cm above the mattress. Depending on whether you have problems of cervical lordosis, or whether you use the pillow prophylactically, it should be turned to the right side. People with cervical hypolordosis are recommended to use the side with the cervical-thoracic profile resembling a tongue. Other people should put their head on the other side of the pillow. This position of the head guarantees ideal alignment of the neck section of the spine. This position provides comfort and relaxation without causing pain in the neck.

Side position

Depending on your height and weight and the firmness of the mattress, the distance between your head and the mattress should be between 9 cm and 11 cm. **The contour of the face** (forehead, nose, mouth, and chin) should protrude slightly beyond the side edge of the pillow. This positioning of the head ensures its proper resting and prevents the head from sliding into the central recess of the pillow. It also provides ideal conditions for breathing smoothly while sleeping.

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


TECHNOLOGY

OVERNIGHT SKIN CARE

Skin+ is a breakthrough in skin care - an innovative fabric treatment based on patent-protected probiotic technology.

Skin+ is a textile material enhanced with carotenoid probiotics with their daily dosage comparable to the effect produced by everyday application of a carotenoid cream.

The Skin+ fabric contains nanocapsules with the proprietary probiotic formula of 'Bacillus indicus Pd01'. As a result of the friction between the body and the fabric, the capsules open up, the probiotics multiply and, as a consequence of the inherent feature of their metabolism, produce carotenoids, which float in the air and get attracted by our skin.

Carotenoids stimulate fibroblasts to synthesize collagen and elastin, resulting in smoother skin wrinkles.

Tested for harmful substances according to Oeko-Tex® Standard 100 Institute

Fabrics of the highest quality,
ECO TEX STANDARD certified.

Foam highest quality,
ISO 9001: 2009 certified.

www.onsensleeping.com

info@onsensleeping.com

+48 509 940 633

ONSEN Sp. z o.o. Marii Skłodowskiej-Curie 71 G, 87-100 Toruń, Poland

The pillow has been entirely designed and manufactured in Poland with the use of the best European materials.